


Bovine Lameness

Foot and Claw Surgery

Instruments
















Gigli Wire



Scalpel



Suture Material and Needle

	Outer Diameter (mm)	Gauge Needle	Hub Color
	0.23	32	Orange
	0.25	31	Violet
	0.30	30	Yellow
	0.33	29	Red
	0.36	28	Blue Green
	0.40	27	Medium Grey
	0.45	26	Brown
	0.50	25	Orange
	0.60	23	Deep Blue
	0.70	22	Black
	0.80	21	Deep Green
	0.90	20	Yellow
	1.10	19	Cream
	1.20	18	Pink
	1.60	16	White



Usually use an 18G to suture skin of cattle

Haemostat

